Vidrio Industrial POBEL, S.A.

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1.1. GENERAL INFORMATION

- The package must be handled carefully. When you unpackage it, you should check that what it contains inside is what is indicated in the Packing List section. If any component is missing or otherwise damaged, contact the dealer immediately.
- Please, read the instructions before installing the apparatus. The instruction manual will be an inseparable part of the product; For this reason, keep it in good condition for future reference.
- The modification of the components of the product transgresses the directive of use 2009/104 / CE; and therefore: the manufacturer is not responsible for the damages that this brings as damage.
- The product must be installed on a horizontal surface and never in places where there are vapors of contaminating or flammable substances.
- Please, dispose of at least the first distillation before using the equipment.
- Conservative packaging during the warranty period.
- Warning and Observation are indicated in you can see below:

Warning: Must be strictly followed to avoid physical damage.

Note: Contains suggestions and restrictions for the product.



1.2 PACKING LIST

Description	Qty.
Metallic water still	1
Feed water tuve	1
Filter	1
Silicone tube 12x18	2
Manual	1

1.3 PRODUCT SPECIFICATONS

Specifications	10DESA0041	10DESA1041	10DESA0081	10DESA0120
Distillation capacity (L/H)	4	4	8	12
Cooling water (L/MIN)	1	1	1,4	1,6
Conductivity 20°C (μS/cm)	2	2	2	2
Wattage (W)	3000	3000	6000	7800
Voltage (V)	230	110	230	380
Intensity (A)	13	27	26	11,3 / phase
Weight (Kg)	8		12	13
Dimensions, height x width x depth (mm)	435 x 370 x 220		635 x 370 x 260	670 x 370 x 260
Safety	Thermostat, level switch			

1.3.1. WATER DISTILLER. FRONT PANEL

1.4. INSTALLATION

1. Connect the equipment to a 230 V power outlet using the power supply cable (110 V for 10DESA1041 and a three-phase connection for 10DESA0120)
2. Connect the gray hose to the water supply socket. The other end of the hose, with the filter previously installed, will be connected to a $\frac{3}{4}$ "tap.
3. The silicone tube must be connected to the outlet for distilled water and the other end to a collecting container whose level is below the outlet by about 25 cm.

1.5. START-UP

1. Check that the installation has been carried out properly.
2. Open the water inlet tap previously connected to a ¼ ”outlet, feed water with a flow of 60 L / H for 10DESA0041 and 10DESA1041. The flow should be 80 L / H for the 10DESA0081 stills; and 100 L / H for 10DESA0120.
3. Turn switch S to the ON position. The green light will come on and after about 10 minutes the distillation process will begin.

Warning: Always disconnect the power supply before carrying out the cleaning operation

1.6. MAINTENANCE

1. With the switch in the OFF position; and the power cable disconnected from the mains remove the refrigerant.
2. Pour in hydrochloric acid, strong water, or hot water until it comes out of the drain. Leave on for 15 minutes.
3. Pour all liquid and impurities into a drainpipe.
4. Rinse with plenty of water.
5. Repeat steps 2, 3 and 4 if necessary.
6. Carry out a distillation that must be discarded.

1.7. SPARE PARTS

Description	Symbol	Cat. Number	Qty.
Contactora, 230 V (10DESA0041 y 10DESA0081)	C1	YE.DCONT00011	1
Contactora, 230 V (10DESA0120)	C1	YE.DCONT00010	1
Contactora, 110 V (10DESA1041)	C1	YE.DCONT10001	1
Termostato	T	YE.DTERM00001	1
Calentador, 110 V (10DESA1041)	R	YE.DRDES00110	1
Calentador, 230 V (10DESA0041 - 10DESA0081)	R - R1,R2	YE.DRDES00220	1 - 2
Calentador, 230 V (10DESA0120)	R1, R2, R3	0100721552	3
Silicona gasket (10DESA0081 y 10DESA0120)	---	YP.DJUSD00075	1
Silicona gasket (10DESA0041 y 10DESA1041)	---	YP.DJUSD00035	1
Plastic switch level	N	YE.DSONI00001	1



1.8. TROUBLE SHOOTING



Problem	Action
The power cord is connected to the mains, the electrical switch S is ON and the green light does not turn on.	Measure and check that the voltage in the network is 230 V. (110 V for 10DESA1041 and 380 v for 10DESA0120) Check the status of the switch.
Water supply is adequate and distillation does not happen	Following this order: Check level switch N. The circuit closes upwards and activates contactor C1. Replace if necessary. There is current supply to the input of contactor C1 and it does not work: replace it. There is current supply to the resistance: R, R1, R2, R3 and it does not heat. Replace the resistance.
In the event of overheating that activates the safety mechanism.	Press the button RESET.


Note:

It will not be allowed to repair any water still that is not properly cleaned.

1.9. ELECTRIC DIAGRAM

1.10. WARRANTY

1. The product covered by this Certificate is protected for two (2) YEARS after the date of the purchase delivery note of the product.
2. During this period of time, in case of manufacturing fault of the product, the user has the right to its servicing, its replacement (except if one of this options was impossible or disproportionate), to have a reduction in the price and to the resolution of the contract (when the purchaser could not require the servicing or the replacement). The servicing will be free of charge for the user. All the servicing and replacements will be carried out in a reasonable period of time.
3. The breakdowns, flaws or caused damages in the product will be excluded as consequence of: abuses, improper use, exposure to adverse weather conditions, spilled batteries, blows, falls, bad or inadequate installation, power failure, waste by normal use, servicing carried out by different people indicated by us, as well as modification of this Certificate and / or another external causes to the product.
4. The glass parts of our distillers are excluded of this warranty.
5. This warranty will be automatically cancelled in case of wrong handling of the product or there is a breach of the manual or instructions regarding the use, assemble and operation.
6. The warranty does not affect to the rights of the customer according to the law.
7. In case the user needs technical help regarding the guaranteed product, the user will contact to the distributor.
8. With the exception of choice of law rules under Spanish private international law and the provisions of the UN Convention on Contracts for the International Sale of Goods (“UN-CISG”), the substantive law of Spain shall apply to all legal relationships and transactions established by this purchase agreement.
9. Please, should give copy of the delivery note or purchase invoice with the Warranty.