**AGITADOR MAGNÉTICO DIGITAL CON CALEFACIÓN**  
*DIGITAL MAGNETIC STIRRER WITH HEATING SYSTEM*

**10AGIT019T**

**MANUAL DE INSTRUCCIONES**

*INSTRUCTIONS MANUAL*

Vidrio Industrial POBEL SA  
C/ Luis I, nº 60 – Almacén 12 – Nave 4 AB1  
Pol. Ind. Vallecas IV – 28031 Madrid (SPAIN)  
T. +34 91 380 33 18 – F. +34 91 380 32 16

<http://www.pobel.com>

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## 1. GENERAL INFORMATION

Pobel heating magnetic stirrer is a laboratory instrument used for liquid mixing. It is mainly used for stirring or simultaneously heating and stirring low-viscosity liquid or solid-liquid mixtures. The shell is aluminum alloy die-casting, beautiful appearance, and the structure can prevent liquid from overflowing into the machine; brushless DC motor and strong magnetic magnetic steel drive, aluminum heating plate, high heating power, fast heating; and it equipped with high temperature protection device, which is widely used Chemical, medical, scientific research, teaching and other fields.

### 1.1. WARRANTY

You have purchased a faithful instrument. This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of twenty-four (24) months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation. For claim under the warranty, please, contact with us. ou may also send the instrument direct to our facilities or we send you the spare parts to help you resolve the problem in next order, enclosing the invoice copy and by giving reasons for the claim. You would be solely liable for freight costs

## 2. LIST OF ITEMS

List of materials contained in the packaging and quantity:

Description	Quantity
Magnetic stirrer	1
Power cord	1
Temperature probe	1
Support for probe	1
Magnetic stirring bar	1
Spare fuse	1
User manual	1

List of items

## 3. TECHNICAL CHARACTERISTICS

Reference	10AGIT019T
Size of the heating plate [cm]	19x19
Material of the heating plate	Ceramic
Max. temperature of the heating plate [°C]	350 ± 10%
Maximum volume [L]	5
Speed in r.p.m. [rpm]	200-2000
Voltage [V]	220 V - 50/60 Hz
Heating power [W]	600

Technical characteristics

## 4. PANEL OPERATION

**TEMP:** set or view the temperature setting value and its internal parameters.  
**SPEED:** set or view the speed setting value and its internal parameters.  
**TIME:** set the running time. After the timing is over, the time setting operation can restart the temperature control operation.

**▲/SHIFT:** In the non-setting state, long press this key for 3 seconds to switch between 2 channels of temperature and save it automatically; in the setting state, click this key to increase the set value, and long press this button to increase the set value continuously.

**▼/Auto-tuning:** In the non-setting state, press this key for 6 seconds to enter the temperature auto-tuning selection state; in the setting state, click this key to decrease the setting value, long press this key to decrease the setting value continuously.

**SWITCH:** on-off key of the stirring (speed).

## 5. OPERATION GUIDE

① Place the machine steadily, turn on the power, and turn on the switch on the side of the machine.

### ② Setting of temperature, speed and time

In the non-setting state, press the **【TEMP】** key to enter the temperature setting mode, the temperature display window value flashes, press the **【▲】** key and **【▼】** key to set the desired value, and then press **【TEMP】** to set the value. Automatically save and return to normal display state.

The setting method of speed and time is the same as above.

### ③ Timing function

Timing adopts countdown mode, when timing, you can choose to run timing mode or constant temperature timing mode; timing time can choose hours or minutes;

When the running time is set to "0", it means that the timing function is not turned on and the machine is in continuous operation mode;

When the running time setting value > 0, it means that the timing function is turned on; if the running timing mode is selected.

### ④ Agitating speed setting

After power on, the speed window displays "OFF", press **【SWITCH】** to enter the stirring speed setting mode.

### ⑤ Temperature control switching function

In the non-set state, long press the 【▲】 key for 3 seconds to switch between the liquid temperature (A mode) and the heating plate surface temperature (B) mode. When it is switched to the liquid temperature control mode, "A" lights up; when switching to the heating plate control mode, "B" lights up.

### ⑥ Over temperature alarm

When the temperature deviation and over-temperature alarm occurs, the buzzer will beep bit by bit, the "ALM" character will light up, and the machine will stop heating.

### ⑦ Abnormal temperature measurement alarm

If "E-X" is displayed in the upper row of the display window, it means that the temperature sensor is faulty or the temperature exceeds the measuring range or the machine is faulty. The machine will automatically stop heating, the buzzer will beep continuously, and the "ALM" character will light up. Please check the temperature sensor and wiring carefully.

#### Note:

E-1: Indicates the failure of the liquid temperature (A mode), that is, the temperature probe failure.

E-2: Indicates that the heating plate surface (B mode) is faulty, that is, the machine is faulty.

E-3: Indicates that the ambient temperature is faulty.

⑧ Press any key to silence the buzzer.

## 6. CALIBRATION

When the temperature control effect is not ideal, the system can be self-tuning. During the auto-tuning process, the temperature will have a large overshoot, and the user is in the system. Please fully consider this factor before auto-tuning.

In the non-setting state, long press the 【▼】 key for 6 seconds to enter the system auto-tuning selection state, and the temperature window will display the auto-tuning prompt. Indicator "AT", the speed window displays the value, you can click the 【▲】 and 【▼】 keys to choose to display "0" or "1".

When "1", click the 【TIME】 key, the controller enters the system auto-tuning state, the [AT] character flashes, after the auto-tuning is completed, [AT] characters stop flashing, the controller will get a better set of system PID parameters, and the parameter values will be saved automatically. Auto-tuned in the system.

During the process, long press 【▼】 key for 6 seconds to stop the auto-tuning program. The "Temperature" key is invalid during system auto-tuning.

## 7. TECHNICAL PROBLEMS

Problem	Action
The machine does not turn on.	Check that it is correctly plugged into the power supply. Check if the fuse is either broken or loose.
An error arises when it is turned on.	Turn off the device, turn it on again and restore it to its default values.
The temperature cannot reach the value adjusted in the settings.	Check in the settings if the established safety temperature is too low.
The stirring speed cannot reach the value adjusted in the settings.	An excessive viscosity of the media may cause a decrease in the speed of the motor.
Heating does not occur after having selected any value of temperature, or there is no stirring after having adjusted it to a specific value.	Check if the control panel has any damage that may have arose during transport.

### Technical problems

If you cannot find a solution to any of the mentioned problems or others, please restore the device to its default values. In necessary, take the unit to the technical service or contact with the manufacturer.

## 8. MAINTENANCE AND CLEANING

- Proper maintenance can keep devices in a good state and can lengthen their lifetime.
- Be careful not to spray the cleanser into the instrument while cleaning.
- Unplug the device from the energy supply while cleaning.
- Use the recommended solutions described in the table to clean the dirt in each case

Type of Dirt	Cleansing Solution
Ink	Isopropyl alcohol (common alcohol)
Construction Materials	Isopropyl alcohol or water with detergent
Cosmetics	Isopropyl alcohol or water with detergent
Food	Water with detergent
Gasoline	Water with detergent

- Wear protective gloves during the cleaning of the instrument.
- Before using other methods for cleansing or decontamination, please contact with the manufacturer to ascertain that this will not cause damage to the device.
- The enamel coating of the heating plate allows its protection and endows it with resistance to acids and bases. However, it also makes it more susceptible to great variations in temperature and to strong impacts, which may provoke the formation of cracks or the detachment of the coating.
- The instrument must be cleaned and put into its initial packaging before sending it to the technical service for repair.
- Use the device in a dry clean room, with a stable temperature.

## 9. SAFETY INDICATION

When working, wear the personal protection to avoid the risk from:

- splashing and evaporation of liquids.
- release of toxic or combustible gases.

Set up the instrument in a spacious area on a stable, clean, non-slip, dry and fireproof surface, do not operate the instrument in explosive atmospheres, with hazardous substances or under water.

Gradually increase the speed, reduce the speed if:

- the stirring bar breakaway because of too high speed.
- the instrument is not running smoothly, or container moves on the base plate.

Temperature must always be set to at least 25 °C lower than the fir point of the media used.

Beware of hazards due to:

- flammable material or media with a low boiling temperature.
- overfilling of media.
- unsafe container.

Processes with pathogenic materials shall be carried out only in closed vessels.

If the stirrer bar is made of PTFE, please note:

- elemental fluorine, three fluoride and alkali metals will corrode the PTFE and halogen alkanes make it expansion at room temperature.
- Molten alkali, alkaline earth metals or their solution, as well as the power in second and third ethnic of the Periodic Table of elements will have chemical reaction with PTFE when temperature reaches 300 – 400 °C.

Check the instrument and accessories before using for damage each time you see them.

- Do not use damaged components. Safe operation is only guaranteed with the accessories described in the “Accessories” chapter.
- Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accesories.

Ensure that the external temperature sensor is inserted in the media to a depth of at least 20 mm.

When using metal vessels, do not place the temperature sensors on the bottom of the vessel.

Placing sensors on the vessel bottom can cause excessively high temperature to be measured especially in media which have poor conductivity. The tip of the measuring sensor must be at least 5 mm from the vessel bottom, a distance of 10 mm is ideal.

The instrument can only be disconnected from the main power supply by pulling out the mains plug or the connector plug.

The voltage stated on the label must correspond to the main power supply.

Make sure that the mains power supply cable does not touch the heating base plate. Do not cover the device. Keep away from high magnetic field. Observe the minimum distances between the devices, between the device and the wall and above the assembly (min. 100 mm).

# POBEL

Vidrio Industrial POBEL SA  
C/ Luis I, nº 60 – Almacén 12 – Nave 4 AB1  
Pol. Ind. Vallecas IV – 28031 Madrid (SPAIN)  
T. +34 91 380 33 18 – F. +34 91 380 32 16  
<http://www.pobel.com>