

Información general

Nombre de producto: CLIPS PINZA PARA ESMERILADO

Fecha Rev.: Dic 04, 2024

Detalles:

Clips de poliacetal (POM) diseñados para juntas esmeriladas, capaces de resistir temperaturas de hasta aproximadamente 150 °C, con excelente resistencia química.

Son fáciles, seguros y rápidos de instalar y retirar.

Cada tamaño se presenta en un color distinto para simplificar su identificación.

No dañan el vidrio.

Mantienen su forma sin doblarse.

Resistentes a la corrosión.

Fotografía



Datos técnicos

Ref.	Para esmerilado NS	Color	Unidades (Pack)
10CSEP1019	10/19	Verde	10
10CSEP1221	12/21	Violeta	10
10CSEP1423	14/23	Amarillo	10
10CSEP1926	19/26	Azul	10
10CSEP2429	24/29	Verde	10
10CSEP2932	29/32	Rojo	10
10CSEP3435	34/35	Naranja	10
10CSEP4540	45/40	Marrón	10



General Information

Product name: CONICAL HEAD CLIP

Rev. Date: Dec 04, 2024

Details:

Ground joint clamps made of polyacetal (POM), heat-resistant up to approx. 150°C, very good chemical resistance.

Easy, safe and quick to fit and remove.

Each size in a different color for easy selection.

No glass damage.

No bending, constant shape.

No corrosion.

Picture



Technical Information

Ref.	Joint	Colour	Units (Pack)
10CSEP1019	10/19	Verde	10
10CSEP1221	12/21	Violet	10
10CSEP1423	14/23	Amarillo	10
10CSEP1926	19/26	Blue	10
10CSEP2429	24/29	Verde	10
10CSEP2932	29/32	Red	10
10CSEP3435	34/35	Naranja	10
10CSEP4540	45/40	Brown	10

