

Información general

Nombre de producto: MATRACES FONDO REDONDO 2 BOCAS PARALELAS

Fecha Rev.: Dic 12, 2024

Detalles:

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Esmerilado P	Esmerilado O
0100104041	100	14/23	14/23
0100104042	100	14/23	19/26
0100104045	100	14/23	29/32
0100104052	250	14/23	19/26
0100104054	250	14/23	29/32
0100104055	250	29/32	29/32
0100104064	500	14/23	29/32
0100104065	500	29/32	29/32
0100104084	1000	14/23	29/32
0100104085	1000	29/32	29/32
0100104094	2000	14/23	29/32
0100104095	2000	29/32	29/32




General Information

Product name: ROUND BOTTOM FLASKS, TWO PARALLEL SIDE NECKS

Rev. Date: Dec 12, 2024

Details:

Picture



Technical Information

Ref.	Capacity (ML)	Joint P	Joint O
0100104041	100	14/23	14/23
0100104042	100	14/23	19/26
0100104045	100	14/23	29/32
0100104052	250	14/23	19/26
0100104054	250	14/23	29/32
0100104055	250	29/32	29/32
0100104064	500	14/23	29/32
0100104065	500	29/32	29/32
0100104084	1000	14/23	29/32
0100104085	1000	29/32	29/32
0100104094	2000	14/23	29/32
0100104095	2000	29/32	29/32

