

Información general


Nombre de producto: MATRACES FONDO REDONDO 3 BOCAS

Fecha Rev.: Dic 10, 2024

Detalles:

Dos bocas paralelas y una inclinada.

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Esmerilado P	Esmerilado N	Esmerilado O
0100107052	250	19/26	14/23	19/26
0100107064	500	19/26	14/23	29/32
0100107084	1000	29/32	29/32	29/32
0100107094	2000	29/32	29/32	29/32
0100107107	3000	29/32	29/32	45/40
0100107117	4000	29/32	29/32	45/40
0100107127	5000	29/32	29/32	45/40
0100107137	6000	29/32	29/32	45/40
0100107147	10000	29/32	29/32	45/40



General Information


Product name: MATRACES FONDO REDONDO 3 BOCAS

Rev. Date: Dec 10, 2024

Details:

Two parallel necks and one inclined neck.

Picture



Technical Information

Ref.	Capacity (ML)	Joint P	Joint N	Joint O
0100107052	250	19/26	14/23	19/26
0100107064	500	19/26	14/23	29/32
0100107084	1000	29/32	29/32	29/32
0100107094	2000	29/32	29/32	29/32
0100107107	3000	29/32	29/32	45/40
0100107117	4000	29/32	29/32	45/40
0100107127	5000	29/32	29/32	45/40
0100107137	6000	29/32	29/32	45/40
0100107147	10000	29/32	29/32	45/40

