

Información general

Nombre de producto: EQUIPO PARA DETERMINACIÓN DE TRIQUINA

Fecha Rev.: Dec 16, 2024

Detalles:

Equipo para la determinación de triquina.

Este equipo está formado por los siguientes componentes:

Ref.	Descripción
10EMTR2000	EMBUDO DECANTACIÓN 2 L
10EMTR0040	EMBUDO CON LLAVE 40 ML
10TAMI0180	TAMIZ ACERO INOX. 100X40 MM 180 UM
10SOPL0900	SOPORTE PIE PLATO 165X300 VARILLA 900 MM
10ACCN0100	ARO CON NUEZ 100X190 MM
10NUDO0110	NUEZ DOBLE P/VARILLA 13 MM
0200913816	VASO PRECIPITADO FORMA BAJA 3000 ML
10PALG0025	PROBETA 25 ML CLASE A
10ASPS0125	ASPIRADOR PIPETAS 0-25 ML
10ESPA0010	ESPÁTULA PP 2X235 MM
10VAMP0608	VARILLA MAGNÉTICA 60X10 MM
10PLRO0010	PLACAS RODAC (25 UNI.)
0200493155	PIPETA GRADUADA 25 ML CLASE B
10PBCR0100	PINZA PARA BURETA Y REFRIGERANTE, SIN NUEZ
10AGIT019T	AGITADOR MAGNÉTICO CON Sonda de TEMP.



Fotografía



Datos técnicos

Ref.	Descripción
10EQTR0001	Equipo completo

el m o l

General Information

Product name: EQUIPMENT FOR TRICHINELLA DIGESTION

Rev. Date: Dec 16, 2024

Details:

Equipment for trichinella digestion.

This kit is made up of the following components:

Ref.	Description
10EMTR2000	SEPARATING FUNNEL 2 L
10EMTR0040	FUNNEL 40 ML
10TAMI0180	TEST SIEVE 100X40 MM. 180 UM
10SOPL0900	RETORT STAND BASE, METALLIC 150 x 300 MM
10ACCN0100	RETORT RINGS W/BOSSHEAD,CLOSED 100X190 MM.
10NUDO0110	BOSSHEAD METALLIC FOR ROD 13MM
0200913816	BEAKER LOW FORM GRADUATED 3000 ML
10PALG0025	GLASS MEASURING CYLINDER 25 ML. CLASS A
10ASPS0125	PIPETTE PUMP 0-25 ML.
10ESPA0010	SPATULA PP 2X235 MM
10VAMP0608	MAGNETIC STIRRING BAR POLYGON 60X10 MM
10PLRO0010	RODAC PLATES (25 PCS)
0200493155	MEASURING PIPETTE GLASS 25 ML. CLASS B
10PBCR0100	CLAMPS FOR BURETTE AND CONDENSERS
10AGIT019T	HOTPLATE STIRRER W/TEMPERATURE SENSOR



Picture



Technical Information

Ref.	Description
10EQTR0001	Complete set