

Información general

Nombre de producto: EMBUDO DE DECANTACIÓN CÓNICO, GRADUADO, LLAVE DE VIDRIO

Fecha Rev.: Feb 07, 2025

Detalles:

Fabricado en vidrio borosilicato 3.3, ofreciendo alta resistencia térmica y química. Cuenta con una junta NS y un tapón de PE, además de un tapón macizo de vidrio para un sellado seguro. Incluye un dispositivo de retención de rosca para mayor estabilidad y una escala graduada que indica el volumen aproximado en mililitros, facilitando la medición de líquidos.

Fotografía



Datos técnicos

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Subdivisiones (ML)
0100239444	19/26	Vidrio	100	5
0100239445	29/32	Vidrio	250	10
0100239446	29/32	Vidrio	500	20
0100239448	29/32	Vidrio	1000	50
0100239450	29/32	Vidrio	2000	100

General Information

Product name: SEPARATING FUNNEL, WITH SCALE, GLASS STOPCOCK

Rev. Date: Feb 07, 2025

Details:

Made of borosilicate glass 3.3, offering high thermal and chemical resistance. It features an NS joint and a PE stopper, as well as a solid glass stopper for a secure seal. It includes a screwthread retaining device for added stability and a graduated scale indicating the approximate volume in milliliters, facilitating liquid measurement.

Picture



Technical Information

Ref.	Stopper	Stopcock	Capacity (ML)	Subdivisiones (ML)
0100239444	19/26	Glass	100 ML	5
0100239445	29/32	Glass	250 ML	10
0100239446	29/32	Glass	500	20
0100239448	29/32	Glass	1000 ML	50
0100239450	29/32	Glass	2000 ML	100

