

Información general

Nombre de producto: EMBUDO DE DECANTACIÓN CÓNICO, GRADUADO, LLAVE PTFE CÓNICA

Fecha Rev.: Feb 07, 2025

Detalles:

Fabricado en vidrio borosilicato 3.3, lo que garantiza alta resistencia térmica y química. Estos embudos cuentan con un tapón de PTFE para un sellado seguro, una escala graduada que indica el volumen aproximado en mililitros, y una junta NS con tapón de PE. Además, incluyen un dispositivo de retención de rosca para mayor estabilidad, asegurando un uso cómodo y preciso en aplicaciones de laboratorio.

Fotografía



Datos técnicos

Ref.	Tapón	Tipo de llave	Capacidad (ML)	Subdivisiones (ML)
0100239544	19/26	PTFE cónica	100	5
0100239545	29/32	PTFE cónica	250	10
0100239546	29/32	PTFE cónica	500	20
0100239548	29/32	PTFE cónica	1000	50
0100239550	29/32	PTFE cónica	2000	100



General Information

Product name: SEPARATING FUNNEL, WITH SCALE, PTFE CONICAL

Rev. Date: Jan 31, 2025

Details:

Made of borosilicate glass 3.3, ensuring high thermal and chemical resistance. These funnels feature a PTFE stopper for a secure seal, a graduated scale indicating the approximate volume in milliliters, and an NS joint with a PE stopper. Additionally, they include a screwthread retaining device for added stability, ensuring comfortable and precise use in laboratory applications.

Picture



Technical Information

Ref.	Stopper	Stopcock	Capacity (ML)	Subdivisiones (ML)
0100239544	19/26	PTFE conical	100 ML	5
0100239545	29/32	PTFE conical	250 ML	10
0100239546	29/32	PTFE conical	500	20
0100239548	29/32	PTFE conical	1000 ML	50
0100239550	29/32	PTFE conical	2000 ML	100

