

Información general

Nombre de producto: TUBOS DE CENTRÍFUGA ASTM, FORMA PERA, GRADUADOS

Fecha Rev.: Feb 12, 2025

Detalles:

Tubo de centrifuga ASTM, forma pera, graduado, con fondo cilíndrico, de acuerdo a ASTM D96. Capacidad de 100 mL, puede ser usado hasta RCF 700 (fuerza centrífuga relativa). Fabricado en vidrio de borosilicato 3.3, utilizado para la determinación de agua y sedimentos en petróleos crudos según el método de ensayo ASTM D96. Presenta elevada resistencia química y térmica.

Fotografía



Datos técnicos

| Ref.       | Designación ASTM | Capacidad (ML) |
|------------|------------------|----------------|
| 0100497200 | D96              | 100            |

## General Information

**Product name:** ASTM CENTRIFUGE TUBES, PEAR SHAPE, GRADUATED

**Rev. Date:** Feb 12, 2025

**Details:**

ASTM centrifuge tubes, pear shape, graduated, with a cylindrical bottom, in accordance with ASTM D96. Capacity of 100 mL, can be used up to RCF 700 (relative centrifugal force). Made of borosilicate glass 3.3, used for the determination of water and sediments in crude oils according to the ASTM D96 test method. It has high chemical and thermal resistance.

## Picture



## Technical Information

| Ref.       | Number ASTM | Capacity (ML) |
|------------|-------------|---------------|
| 0100497200 | D96         | 100 ML        |

