

Información general

Nombre de producto: TUBOS DE CENTRÍFUGA, FONDO REDONDO, TAPÓN DE BAQUELITA

Fecha Rev.: Feb 11, 2025

Detalles:

Fabricado en vidrio de borosilicato 3.3 de pared gruesa, este tubo de centrífuga está diseñado para resistir altas velocidades y garantizar máxima durabilidad. Su composición le otorga una excelente resistencia a productos químicos agresivos y cambios bruscos de temperatura. Es autoclavable y puede esterilizarse por calor seco y vapor hasta 180°C.

Cuenta con un diseño cilíndrico de fondo redondo y un tapón de baquelita a rosca con interior de PTFE. El sistema de sellado incluye una rosca SVL y un disco de etileno propileno recubierto de PTFE, asegurando una protección óptima. Soporta fuerzas de hasta 3500 ×g, ofreciendo un rendimiento fiable en aplicaciones de laboratorio exigentes.

Fotografía



Datos técnicos

Ref.	Ø Exterior (MM)	Longitud (mm)	Rosca	Capacidad (ML)
0100159200	16	110	SVL-13	110
0100159202	24	110	SVL-18	110
0100159204	34	110	SVL-25	110
0100159206	40	130	SVL-30	130
0100159208	44	130	SVL-30	130
0100159210	56	165	SVL-42	165



General Information

Product name: CENTRIFUGE TUBES, ROUND BOTTOM, BAKELITE SCREW CAP

Rev. Date: Feb 11, 2025

Details:

Made of thick-walled borosilicate glass 3.3, these centrifuge tubes are designed to withstand high speeds and ensure maximum durability. Their composition provides excellent resistance to aggressive chemicals and sudden temperature changes. They are autoclavable and can be sterilized by dry heat or steam up to 180°C.

Featuring a cylindrical design with a round bottom, they come with a bakelite screw cap lined with PTFE. The sealing system includes an SVL screw thread and an ethylene-propylene disc coated with PTFE, ensuring optimal protection. Withstanding forces up to 3500 xg, they offer reliable performance for demanding laboratory applications.

Picture



Technical Information

Ref.	Ø Exterior (MM)	Length (mm)	Thread	Capacity (ML)
0100159200	16	110	SVL-13	110
0100159202	24	110	SVL-18	110
0100159204	34	110	SVL-25	110
0100159206	40	130	SVL-30	130
0100159208	44	130	SVL-30	130
0100159210	56	165	SVL-42	165

