

Información general

Nombre de producto: TUBOS DE CENTRÍFUGA CON ROSCA, FONDO REDONDO, TAPÓN PBT

Fecha Rev.: Feb 12, 2025

Detalles:

Los tubos de centrífuga con fondo redondo y tapones de rosca están diseñados para soportar altas velocidades de centrifugado. Fabricados con paredes más gruesas, ofrecen una mayor resistencia y durabilidad. Cuentan con tapones de PBT y una junta de silicona/PTFE que garantiza un cierre seguro y hermético.

Son completamente autoclavables, lo que permite su esterilización sin comprometer su integridad. Están fabricados en vidrio que cumple con los requisitos de la USP Tipo I y la norma ASTM E438, Tipo I, Clase A. Su fondo es redondo y el sistema de cierre es mediante un tapón de rosca. El material utilizado es borosilicato 3.3, reconocido por su resistencia química y térmica.

Fotografía



Datos técnicos

Ref.	Ø Exterior (MM)	Longitud (mm)	Rosca	Capacidad (ML)
0100159300	16	100	GL-18	10
0100159302	24	110	GL-25	22
0100159304	34	110	GL-25	45
0100159306	40	130	GL-32	75
0100159308	44	130	GL-32	95



General Information

Product name: CENTRIFUGE TUBES WITH THREAD, ROUND BOTTOM, PBT CAP.

Rev. Date: Feb 2, 2025

Details:

Round bottom centrifuge tubes with screw caps are designed to withstand high centrifugation speeds. Manufactured with thicker walls, they offer greater strength and durability. They feature PBT caps with a silicone/PTFE seal, ensuring a secure and airtight closure.

These tubes are fully autoclavable, allowing for sterilization without compromising their integrity. They are made of glass that meets USP Type I and ASTM E438, Type I, Class A standards. The bottom is round, and the closure system consists of a screw cap. The material used is borosilicate 3.3, known for its chemical and thermal resistance.

Picture



Technical Information

Ref.	Ø Ext. (MM)	Length (mm)	Thread	Capacity (ML)
0100159300	16	100	GL-18	10
0100159302	24	110	GL-25	22
0100159304	34	110	GL-25	45
0100159306	40	130	GL-32	75
0100159308	44	130	GL-32	95

