

Información general

Nombre de producto: MATRACES DE DECANTACIÓN ÁMBAR CON LLAVE DE VIDRIO

Fecha Rev.: Feb 12, 2025

Detalles:

Embudos de decantación según Squibb, de color ámbar, fabricados en vidrio de borosilicato 3.3. Equipados con tapón macizo de vidrio transparente, sin escala. Incorporan encaje normalizado (NS) y tapón de PE, además de un dispositivo de retención de rosca para mayor seguridad.

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Esmerilado Hembra
0100239402AM	25	14/23
0100239403AM	50	14/23
0100239404AM	100	19/26
0100239405AM	250	29/32
0100239406AM	500	29/32
0100239408AM	1000	29/32



General Information

Product name: AMBER DECANTING FLASKS WITH GLASS STOPCOCK

Rev. Date: Feb 12, 2025

Details:

Amber decanting flasks according to Squibb, amber-colored, made of borosilicate glass 3.3. Equipped with a solid clear glass plug, without scale. They feature a standardized NS socket and a PE stopper, along with a screw-thread retaining device for added security.

Picture



Technical Information

Ref.	Capacity (ML)	Socket Joint
0100239402AM	25	14/23
0100239403AM	50	14/23
0100239404AM	100	19/26
0100239405AM	250	29/32
0100239406AM	500	29/32
0100239408AM	1000	29/32

