

Información general


Nombre de producto: SCHILLING BURETTES, CLASS A, COMPLETE SET - AMBER

Fecha Rev.: Jul 11, 2023

Detalles:

- Con enrase automático a cero.
- Circuito cerrado para protección de la contaminación exterior.
- Unidad compacta que hace que no sean necesarios accesorios metálicos (pinzas, etc.).
- Fácil carga de la bureta, simplemente por presión del frasco.

Fotografía



Datos técnicos

Ficha de producto

Ref.	Llave	Capacidad (ML)	Tolerancia (±ML)	Subdivisiones (ML)
0100491478AM	PTFE Aguja	25	± 0.05	0.1
0100491479AM	PTFE Aguja	50	± 0.05	0.1
0100491481AM	Plástico	25	± 0.05	0.1
0100491482AM	Plástico	50	± 0.05	0.1



General Information


Product name: SCHILLING BURETTES, CLASS A, COMPLETE SET - AMBER

Rev. Date: Jul 11, 2023

Details:

- With automatic zeroing.
- Closed circuit to avoid external contamination.
- Compact unit with no need for metallic accesories (clamps, etc.).
- Easy loading from flask pressure.

Picture



Technical Information

Product datasheet

Ref.	Stopcock	Capacity (ML)	Tolerance (±ML)	Subdivisiones (ML)
0100491478AM	PTFE NEEDLE	25 ML	± 0.05	0.1
0100491479AM	PTFE NEEDLE	50 ML	± 0.05	0.1
0100491481AM	Plastic	25 ML	± 0.05	0.1
0100491482AM	Plastic	50 ML	± 0.05	0.1

