

Información general

Nombre de producto: MICROPIPETAS MULTICANAL DE VOLUMEN VARIABLE - HTL MODELO DISCOVERY PRO

Fecha Rev.: Jul 11, 2023

Detalles:

Disponible en 8 y 12 canales.

Calidad acorde a estándares GLP.

Bloqueo del ajuste de volumen.

Código de color.

Contador con 4 dígitos.

Fácil sistema de calibración.

Fotografía



Datos técnicos

Ficha de producto

Ref.	Rango Volumen (μ L)	Color	Canales
154951DP8-10	1-10	Rojo	8
154952DP8-50	5-50	Amarillo	8
154953DP8-200	20-200	Amarillo	8
154954DP8-300	50-300	Verde	8
154955DP12-10	1-10	Rojo	12
154956DP12-50	5-50	Amarillo	12
154957DP12-200	20-200	Amarillo	12
154958DP12-300	50-300	Verde	12



General Information

Product name: HTL MULTICHANNEL ADJUSTABLE VOLUME MICROPIPETTE, MODEL DISCOVERY PRO

Rev. Date: Jul 11, 2023

Details:

Available in 8 or 12 channels.

Quality according to GLP standards.

Volume string lock.

Colour coding.

Four digit counter.

User friendly calibration system.

Picture



Technical Information

Product datasheet

Ref.	Range Volume (μ L)	Colour	Channels
154951DP8-10	1-10	Red	8
154952DP8-50	5-50	Amarillo	8
154953DP8-200	20-200	Amarillo	8
154954DP8-300	50-300	Verde	8
154955DP12-10	1-10	Red	12
154956DP12-50	5-50	Amarillo	12
154957DP12-200	20-200	Amarillo	12
154958DP12-300	50-300	Verde	12

