

Información general

Nombre de producto: MICROPIPETAS MONOCANAL DE VOLUMEN VARIABLE - HTL MODELO DISCOVERY PRO

Fecha Rev.: Jul 11, 2023

Detalles:

Ergonomía, precisión y exactitud.

Correcto agarre, reduciendo el estrés en la muñeca y fatiga, causantes de lesión por esfuerzo repetitivo.

Elimina el riesgo de transferencia de calor al mecanismo interior de la pipeta.

Estándares de calidad GLP.

Bloqueo del ajuste de volumen.

Código de color

Contador de 4 dígitos visible durante el trabajo.

Fácil sistema de calibración.

Fotografía



Datos técnicos

Ref.	Rango Volumen (μ L)	Color
156001DP2	0.2-2.0	Rojo
156002DP10	0,5-10	Rojo
156003DP20	2-20	Amarillo
156004DP100	10-100	Amarillo
156005DP200	20-200	Amarillo
156006DP1000	100-1000	Azul



General Information

Product name: HTL SINGLE CHANNEL ADJUSTABLE VOLUME MICROPIPETTE, MODEL DISCOVERY PRO

Rev. Date: Jul 11, 2023

Details:

Ergonomics, accuracy and precision.

Its ergonomic grip reduces stress on the wrist and fatigue, thus preventing repetitive strain injuries.

Eliminates the risk of heat transfer to the inner pipette mechanism.

Quality according to GLP standards.

Volume string lock.

Colour coding.

Four digit counter remains visible during work.

User friendly calibration system.

Picture



Technical Information

Ref.	Range Volume (μL)	Colour
156001DP2	0.2-2.0	Red
156002DP10	0.5-10	Red
156003DP20	2-20	Amarillo
156004DP100	10-100	Amarillo
156005DP200	20-200	Amarillo
156006DP1000	100-1000	Blue

