

## Información general

**Nombre de producto:** PIPETAS GRADUADAS PARA CULTIVOS

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

Pipetas graduadas para cultivos. Calibradas por vertido (TD, EX).

Con boca de aspiración normalizada para alojamiento de tapón de algodón.

Fabricadas en vidrio soda.

## Fotografía



## Datos técnicos

Ref.	Longitud "L" (MM)	Capacidad (ML)	Subdivisiones (ML)
0200493181	180	1	0.02

## Ficha de producto

Ref.	Longitud &quot;L&quot; (MM)	Capacidad (ML)	Subdivisiones (ML)
0200493182	180	2	0.05
0200493185	180	5	0.1
0200493190	180	10	0.2
0200493195	230	25	0.5



## General Information

**Product name:** GRADUATED PIPETTES FOR TISSUE CULTURE

**Rev. Date:** Jul 11, 2023

**Details:**

Graduated pipette for tissue culture. Calibrated to deliver (TD, EX).

Enlarged end for the use of cotton plugs.

Made from soda lime glass.

## Picture



## Technical Information

Ref.	Length &quot;L&quot; (mm)	Capacity (ML)	Subdivisiones (ML)
0200493181	180	1	0.02

## Product datasheet

Ref.	Length &quot;L&quot; (mm)	Capacity (ML)	Subdivisiones (ML)
0200493182	180	2	0.05
0200493185	180	5 ML	0.1
0200493190	180	10 ML	0.2
0200493195	230	25 ML	0.5

