

Información general

Nombre de producto: ADAPTADOR M/M ACODADO A 90º

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Esmerilado Macho
0100220211	14/23
0100220222	19/26
0100220233	24/29
0100220244	29/32

General Information

Product name: ADAPTER 2 CONES, BENT 90°

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Cone Joint
0100220211	14/23
0100220222	19/26
0100220233	24/29
0100220244	29/32