

Información general

Nombre de producto: ADAPTADOR 2 CUELLOS ACODADO 105°

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Esmerilado Macho Inferior	Esmerilado Hembra Lateral	Esmerilado Hembra
0100222811	14/23	14/23	14/23
0100222822	19/26	19/26	19/26
0100222841	29/32	14/23	29/32
0100222844	29/32	29/32	29/32

General Information

Product name: ADAPTER 2 NECKS, BENT 105°

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Bottom Cone Joint	Side Socket (NS)	Socket Joint
0100222811	14/23	14/23	14/23
0100222822	19/26	19/26	19/26
0100222841	29/32	14/23	29/32
0100222844	29/32	29/32	29/32