

Información general

Nombre de producto: BOLA KJELDAHL INVERTIDA

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Esmerilado Macho	Esmerilado Hembra
0100223111	14/23	14/23
0100223122	19/26	19/26
0100223144	29/32	29/32

General Information

Product name: SPLASH HEAD, ROTATY EVAPORATOR, WITH SINTER

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Cone Joint	Socket Joint
0100223111	14/23	14/23
0100223122	19/26	19/26
0100223144	29/32	29/32