

## Información general

**Nombre de producto:** ADAPTADOR PARA DISTRIBUCIÓN, RECTO

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

Fabricado en vidrio borosilicato 3.3. Rama larga, ideal para la dosificación de líquidos desde los condensadores a los recipientes contenedores.

## Fotografía



## Datos técnicos

Ref.	Esmerilado	Longitud vástago (MM)
0100225101	14/23	150
0100225102	19/26	150

## Ficha de producto

Ref.	Esmerilado	Longitud vástago (MM)
0100225103	24/29	150
0100225104	29/32	150

The image features a large, vertical watermark logo on the right side of the page. The logo consists of the word "ESMERAL" in a bold, stylized, light green font. The letters are thick and blocky, with a slight shadow effect. The 'E' is particularly large and prominent. The logo is oriented vertically, reading from top to bottom.

## General Information

**Product name:** DELIVERY ADAPTER, STRAIGHT

**Rev. Date:** Jul 11, 2023

**Details:**

Produced in borosilicate 3.3 glass. Long stem, useful for distribution of liquid from the condensers into the receiving vessels.

## Picture



## Technical Information

Ref.	Joint	Stem Length (MM)
0100225101	14/23	150
0100225102	19/26	150

## Product datasheet

Ref.	Joint	Stem Length (MM)
0100225103	24/29	150
0100225104	29/32	150

