

Información general

Nombre de producto: SEPARADOR POR GIRO 3 BOCAS

Fecha Rev.: Jul 11, 2023

Detalles:

Para tres matraces colectores de hasta 100 mL.

Fotografía



Datos técnicos

Ref.	Esmerilados Macho	Esmerilado Hembra
0100227311	14/23	14/23
0100227321	14/23	19/26
0100227331	14/23	29/32

General Information

Product name: DISTILLING RECEIVER, 3 NECKS W/OLIVE

Rev. Date: Jul 11, 2023

Details:

For three receiver flasks of up to 100 mL.

Picture



Technical Information

Ref.	Cone Joints	Socket Joint
0100227311	14/23	14/23
0100227321	14/23	19/26
0100227331	14/23	29/32