

Información general

Nombre de producto: MATRACES FONDO REDONDO 3 BOCAS PARALELAS

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Esmerilado P	Esmerilado N	Esmerilado O
0100109001	100	19/26	19/26	24/29
0100109002	100	14/23	14/23	29/32
0100109003	100	19/26	19/26	29/32
0100109004	100	29/32	29/32	29/32

Ficha de producto

Ref.	Capacidad (ML)	Esmerilado P	Esmerilado N	Esmerilado O
0100109005	250	19/26	19/26	24/29
0100109006	250	14/23	14/23	29/32
0100109007	250	19/26	19/26	29/32
0100109008	250	29/32	29/32	29/32
0100109009	500	19/26	19/26	24/29
0100109010	500	14/23	14/23	29/32
0100109011	500	19/26	19/26	29/32
0100109012	500	29/32	29/32	29/32
0100109013	1000	19/26	19/26	24/29
0100109014	1000	14/23	14/23	29/32
0100109015	1000	19/26	19/26	29/32
0100109016	1000	29/32	29/32	29/32
0100109017	1000	29/32	29/32	45/40
0100109018	2000	19/26	19/26	24/29
0100109019	2000	14/23	14/23	29/32
0100109020	2000	19/26	19/26	29/32
0100109021	2000	29/32	29/32	29/32
0100109022	2000	29/32	29/32	45/40
0100109023	3000	14/23	14/23	29/32
0100109024	3000	29/32	29/32	29/32
0100109025	3000	29/32	29/32	45/40
0100109026	4000	14/23	14/23	29/32
0100109027	4000	29/32	29/32	29/32
0100109028	4000	29/32	29/32	45/40
0100109029	5000	29/32	29/32	29/32



Ficha de producto

Ref.	Capacidad (ML)	Esmerilado P	Esmerilado N	Esmerilado O
0100109030	5000	29/32	29/32	45/40
0100109031	6000	14/23	14/23	29/32
0100109032	6000	29/32	29/32	29/32
0100109033	6000	29/32	29/32	45/40
0100109034	10000	29/32	29/32	29/32
0100109035	10000	29/32	29/32	45/40
0100109036	20000	29/32	29/32	45/40



General Information

Product name: MATRACES FONDO REDONDO 3 BOCAS PARALELAS

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Capacity (ML)	Joint P	Joint N	Joint O
0100109001	100 ML	19/26	19/26	24/29
0100109002	100 ML	14/23	14/23	29/32
0100109003	100 ML	19/26	19/26	29/32
0100109004	100 ML	29/32	29/32	29/32

Product datasheet

Ref.	Capacity (ML)	Joint P	Joint N	Joint O
0100109005	250 ML	19/26	19/26	24/29
0100109006	250 ML	14/23	14/23	29/32
0100109007	250 ML	19/26	19/26	29/32
0100109008	250 ML	29/32	29/32	29/32
0100109009	500	19/26	19/26	24/29
0100109010	500	14/23	14/23	29/32
0100109011	500	19/26	19/26	29/32
0100109012	500	29/32	29/32	29/32
0100109013	1000 ML	19/26	19/26	24/29
0100109014	1000 ML	14/23	14/23	29/32
0100109015	1000 ML	19/26	19/26	29/32
0100109016	1000 ML	29/32	29/32	29/32
0100109017	1000 ML	29/32	29/32	45/40
0100109018	2000 ML	19/26	19/26	24/29
0100109019	2000 ML	14/23	14/23	29/32
0100109020	2000 ML	19/26	19/26	29/32
0100109021	2000 ML	29/32	29/32	29/32
0100109022	2000 ML	29/32	29/32	45/40
0100109023	3000 ML	14/23	14/23	29/32
0100109024	3000 ML	29/32	29/32	29/32
0100109025	3000 ML	29/32	29/32	45/40
0100109026	4000	14/23	14/23	29/32
0100109027	4000	29/32	29/32	29/32
0100109028	4000	29/32	29/32	45/40
0100109029	5000	29/32	29/32	29/32



Product datasheet

Ref.	Capacity (ML)	Joint P	Joint N	Joint O
0100109030	5000	29/32	29/32	45/40
0100109031	6000	14/23	14/23	29/32
0100109032	6000	29/32	29/32	29/32
0100109033	6000	29/32	29/32	45/40
0100109034	10000	29/32	29/32	29/32
0100109035	10000	29/32	29/32	45/40
0100109036	20000	29/32	29/32	45/40

