

Información general

Nombre de producto: MATRACES PARA CULTIVO, CUELLO RECTO


Fecha Rev.: Jul 11, 2023

Detalles:

Ideales para su uso en aplicaciones de cultivo celular, estos matraces cónicos de cuello estrecho cuentan con 3 hendiduras deflectoras que proporcionan una mayor agitación del medio de cultivo para mejorar la transferencia de oxígeno o gas cuando se utilizan con agitadores rotatorios o alternativos

- Ideales para su uso en la preparación de cultivos
- Fabricado en vidrio de borosilicato 3.3.
- Presenta 3 hendiduras deflectoras para proporcionar una agitación óptima y consistente
- Apto para ciclos repetidos de autoclave a 121°C

Fotografía



Datos técnicos

Ficha de producto

Ref.	Capacidad (ML)
0100126010	100
0100126025	250
0100126050	500
0100126100	1000



General Information

Product name: ERLLENMEYER FLASKS FOR CULTURE, STRAIGHT NECK


Rev. Date: Jul 11, 2023

Details:

Ideal for use in cell culture applications, these narrow neck conical flasks feature 3 baffle indents that provide greater agitation of the culture media to improve oxygen or gas transfer when used with rotary or reciprocating shakers

- Ideal for use in preparation of cultures
- Manufactured from borosilicate glass 3.3
- Exhibits 3 baffle indents to provide optimal and consistent agitation
- Suitable for repeated autoclave cycles at 121°C

Picture



Technical Information

Product datasheet

Ref.	Capacity (ML)
0100126010	100 ML
0100126025	250 ML
0100126050	500
0100126100	1000 ML

