

Información general

Nombre de producto: MATRACES FORMA CORAZÓN CON 2 BOCAS

Fecha Rev.: Jul 11, 2023

Detalles:

En vidrio borosilicato.

Fotografía



Datos técnicos

Ref.	Capacidad (ML)	Esmerilado P	Esmerilado O
0100133021	25	14/23	14/23
0100133022	25	14/23	19/26
0100133031	50	14/23	14/23

Ficha de producto

Ref.	Capacidad (ML)	Esmerilado P	Esmerilado O
0100133032	50	14/23	19/26
0100133041	100	14/23	14/23
0100133042	100	14/23	19/26
0100133051	250	14/23	14/23
0100133052	250	14/23	19/26



General Information

Product name: PEAR-SHAPED FLASK TWO NECK,

Rev. Date: Jul 11, 2023

Details:

With joint borosilicate glass.

Picture



Technical Information

Ref.	Capacity (ML)	Joint P	Joint O
0100133021	25 ML	14/23	14/23
0100133022	25 ML	14/23	19/26
0100133031	50 ML	14/23	14/23

Product datasheet

Ref.	Capacity (ML)	Joint P	Joint O
0100133032	50 ML	14/23	19/26
0100133041	100 ML	14/23	14/23
0100133042	100 ML	14/23	19/26
0100133051	250 ML	14/23	14/23
0100133052	250 ML	14/23	19/26

