

Información general

Nombre de producto: DESECADORES EN VIDRIO

Fecha Rev.: Jul 11, 2023

Detalles:

Vidrio borosilicato, tapa botón a rosca, sin placa.

Fotografía



Datos técnicos

Ref.	Ø Interior (MM)	Capacidad Aprox. (ML)
0200929210	150	1500
0200929220	200	4000
0200929250	250	6000
0200929260	300	10000

General Information

Product name: GLASS DESICCATORS

Rev. Date: Jul 11, 2023

Details:

Borosilicate glass, button screw cover, without plate.

Picture



Technical Information

Ref.	Ø ID (MM)	Capacity Approx. (ML)
0200929210	150	1500
0200929220	200	4000
0200929250	250	6000
0200929260	300	10000