

Información general

Nombre de producto: REFRIGERANTE DAVIES

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Longitud útil (MM)	Superf. intercambio (CM2)	Esmerilado Hembra/Macho
0100231022	200	250	19/26
0100231034	300	450	29/32

General Information

Product name: DAVIES CONDENSER

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Useful length (MM)	Surface area (CM2)	Socket/Cone Joints
0100231022	200	250	19/26
0100231034	300	450	29/32