


Información general

Nombre de producto: REFRIGERANTE DE SERPENTÍN DIMROTH

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Longitud útil (MM)	Superf. intercambio (CM2)	Esmerilado Hembra/Macho
0100233511	120	140	14/23
0100233522	200	250	19/26
0100233524	250	380	29/32
0100233544	400	600	29/32


General Information

Product name: DIMROTH CONDENSER WITH GLASS OLIVES

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Useful length (MM)	Surface area (CM2)	Socket/Cone Joints
0100233511	120	140	14/23
0100233522	200	250	19/26
0100233524	250	380	29/32
0100233544	400	600	29/32