


Información general

Nombre de producto: REFRIGERANTE DIMROTH PARA REFLUJO Y EXTRACTORES SOXHLET

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Longitud útil (MM)	Esmerilado Macho
0100233614	150	29/32
0100233627	200	45/40
0100233629	200	55/44

General Information

Product name: DIMROTH CONDENSER FOR REFLUX AND SOXHLET EXTRACTION

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Useful length (MM)	Cone Joint
0100233614	150	29/32
0100233627	200	45/40
0100233629	200	55/44