

## Información general

Nombre de producto: REFRIGERANTE DE BOLAS ALLIHN SIN UNIÓN ESMERILADA

Fecha Rev.: Jul 11, 2023

Detalles:

## Fotografía



## Datos técnicos

Ref.	Longitud útil (MM)	Nº Bolas
0100234601	200	4
0100234602	300	6
0100234604	500	10

## General Information

**Product name:** ALLIHN CONDENSER WITHOUT JOINTS

**Rev. Date:** Jul 11, 2023

**Details:**

## Picture



## Technical Information

Ref.	Useful length (MM)	Number of bulbs
0100234601	200	4
0100234602	300	6
0100234604	500	10