

Información general

Nombre de producto: EQUIPO DE FILTRACIÓN CON FRASCO Y SOPORTE FILTRANTE CON OLIVA

Fecha Rev.: Jul 11, 2023

Detalles:

Compuesto por los siguientes productos:

- Cuerpo con placa porosa para anillo corredizo, con oliva de vidrio (Ref. 0100259940) y plástico (Ref. 0100259941)
- Embudo de filtración de 300mL (Ref. 0100260004)
- Pinza de sujeción (Ref. 0100260006)
- Tapón con orificio y rosca GL-45 (Ref. 0200162040)
- Junta de silicona (Ref. 0200162126)
- Frasco ISO graduado de 1000mL

Fotografía



Datos técnicos

Ref.	Oliva	Capacidad (ML)
0100259880	Vidrio	1000
0100259881	Plástico	1000



General Information

Product name: FILTRATION SET WITH BOTTLE AND FILTER HOLDER WITH OLIVE

Rev. Date: Jul 11, 2023

Details:

Includes the following products:

- Body with fritted disc for sliding o-ring. With glass (Ref. 0100259940) and plastic (Ref. 0100259941) olive
- Filter funnel, 300mL (Ref. 0100260004)
- Clamp (Ref. 0100260006)
- Screw cap with aperture, GL-45 thread (Ref. 0200162040)
- Silicone seal (Ref. 0200162126)
- Laboratory bottle, 1000mL

Picture



Technical Information

Ref.	Hose	Capacity (ML)
0100259880	Glass	1000 ML
0100259881	Plastic	1000 ML

