

Información general

Nombre de producto: EQUIPO DE FILTRACIÓN CON KITASATO, TAPÓN Y SOPORTE

Fecha Rev.: Jul 11, 2023

Detalles:

Compuesto por los siguientes productos:

- Embudo con placa filtrante, porosidad nº3 (Ref. 0100260563)
- Cuerpo con placa porosa (Ref. 0100260005)
- Pinza de sujeción (Ref. 0100260006)
- Tapón de silicona con orificio (Ref. 0200259950)
- Matraz kitasato de 1000mL, con oliva de plástico (Ref. 0200913213) y vidrio (Ref. 0200913214)

Fotografía



Datos técnicos

Ficha de producto

Ref.	Oliva	Capacidad (ML)
0100259905	Plástico	1000
0100259906	Vidrio	1000



General Information

Product name: FILTRATION SET WITH FILTERING FLASK, STOPPER AND FILTER HOLDER

Rev. Date: Jul 11, 2023

Details:

Includes the following products:

- Filter funnel with fritted disc, porosity nº3 (Ref. 0100260563)
- Body with fritted disc (Ref. 0100260005)
- Clamp (Ref. 0100260006)
- Silicone stopper with hole (Ref. 0200259950)
- Filter flask, 1000mL, with plastic (Ref. 0200913213) and glass (Ref. 0200913214) olive

Picture



Technical Information

Product datasheet

Ref.	Hose	Capacity (ML)
0100259905	Plastic	1000 ML
0100259906	Glass	1000 ML

