

Información general

Nombre de producto: EMBUDOS DE FILTRACIÓN CON PLACA, FORMA ALEMANA

Fecha Rev.: Jul 11, 2023

Detalles:

Porosidad nº 1 al 4.

Fotografía



Datos técnicos

Ref.	Ø Máximo (MM)	Ø Placa (MM)	Porosidad	Ø Exterior vástago (MM)
0100260751	55	15	1	8
0100260752	55	15	2	8
0100260753	55	15	3	8

Ficha de producto

Ref.	Ø Máximo (MM)	Ø Placa (MM)	Porosidad	Ø Exterior vástago (MM)
0100260754	55	15	4	8
0100260771	70	20	1	8
0100260772	70	20	2	8
0100260773	70	20	3	8
0100260774	70	20	4	8
0100260781	80	25	1	10
0100260782	80	25	2	10
0100260783	80	25	3	10
0100260784	80	25	4	10
0100260791	100	30	1	10
0100260792	100	30	2	10
0100260793	100	30	3	10
0100260794	100	30	4	10



General Information

Product name: FILTER FUNNEL WITH FRITTED DISC, CONICAL SHAPE

Rev. Date: Jul 11, 2023

Details:

Porosity nº 1 to 4.

Picture



Technical Information

Ref.	Maximum \varnothing (MM)	\varnothing Disc (MM)	Porosity	OD Stem (MM)
0100260751	55	15	1	8
0100260752	55	15	2	8
0100260753	55	15	3	8

Product datasheet

Ref.	Maximum \varnothing (MM)	\varnothing Disc (MM)	Porosity	OD Stem (MM)
0100260754	55	15	4	8
0100260771	70	20	1	8
0100260772	70	20	2	8
0100260773	70	20	3	8
0100260774	70	20	4	8
0100260781	80	25	1	10
0100260782	80	25	2	10
0100260783	80	25	3	10
0100260784	80	25	4	10
0100260791	100	30	1	10
0100260792	100	30	2	10
0100260793	100	30	3	10
0100260794	100	30	4	10

