

Información general

Nombre de producto: REACTORES CILÍNDRICOS FONDO REDONDO CON LLAVE

Fecha Rev.: Jul 11, 2023

Detalles:

Cilíndricos; sin graduar. Con llave de descarga.

Fotografía



Datos técnicos

Ref.	Ø Ext. Reactor (mm)	Ø Int. Reactor (mm)	Llave (mm)	Capacidad (ML)	Altura Total (mm)	DN Boca
0100802102	65	60	10	100	185	60
0100802104	65	60	10	250	250	60

Ficha de producto

Ref.	Ø Ext. Reactor (mm)	Ø Int. Reactor (mm)	Llave (mm)	Capacidad (ML)	Altura Total (mm)	DN Boca
0100802106	65	60	10	500	320	60
0100802108	65	60	10	1000	320	60
0100802110	105	99	10	250	260	100
0100802112	105	99	10	500	260	100
0100802114	105	99	10	1000	305	100
0100802116	105	99	10	2000	440	100
0100802118	160	150	10	3000	410	100
0100802120	160	150	10	4000	460	100
0100802121	160	150	10	5000	460	100
0100802122	160	150	10	2000	310	150
0100802124	160	150	10	3000	360	150
0100802126	160	150	10	4000	435	150
0100802128	160	150	10	5000	460	150
0100802130	160	150	10	6000	490	150
0100802132	250	240	20	10000	520	150
0100802134	250	240	20	20000	740	150
0100802136	215	200	10	6000	410	200
0100802138	215	200	20	10000	555	200

General Information

Product name: CYLINDRICAL REACTION VESSELS, ROUND BOTTOM, WITH STOPCOCK

Rev. Date: Jul 11, 2023

Details:

Cylindrical; not graduated. With drain stopcock.

Picture



Technical Information

Ref.	Vessel Outer Ø (mm)	Vessel Inner Ø (mm)	Valve (mm)	Capacity (ML)	Total Height (mm)	Laboratory Flange DN
0100802102	65	60	10	100 ML	185	60
0100802104	65	60	10	250 ML	250	60
0100802106	65	60	10	500	320	60

Product datasheet

Ref.	Vessel Outer Ø (mm)	Vessel Inner Ø (mm)	Valve (mm)	Capacity (ML)	Total Height (mm)	Laboratory Flange DN
0100802108	65	60	10	1000 ML	320	60
0100802110	105	99	10	250 ML	260	100
0100802112	105	99	10	500	260	100
0100802114	105	99	10	1000 ML	305	100
0100802116	105	99	10	2000 ML	440	100
0100802118	160	150	10	3000 ML	410	100
0100802120	160	150	10	4000	460	100
0100802121	160	150	10	5000	460	100
0100802122	160	150	10	2000 ML	310	150
0100802124	160	150	10	3000 ML	360	150
0100802126	160	150	10	4000	435	150
0100802128	160	150	10	5000	460	150
0100802130	160	150	10	6000	490	150
0100802132	250	240	20	10000	520	150
0100802134	250	240	20	20000	740	150
0100802136	215	200	10	6000	410	200
0100802138	215	200	20	10000	555	200