

Información general

Nombre de producto: TAPAS REACTORES 2 BOCAS

Fecha Rev.: Jul 11, 2023

Detalles:

Una boca central vertical y una boca lateral inclinada.

Fotografía



Datos técnicos

Ref.	Altura (mm)	Boca Central	Boca Lateral	DN Boca
0100806030	75	14/23	14/23	60
0100806032	85	29/32	14/23	60
0100806034	95	29/32	14/23	100

Ficha de producto

Ref.	Altura (mm)	Boca Central	Boca Lateral	DN Boca
0100806036	95	29/32	29/32	100
0100806038	110	45/40	29/32	100
0100806040	125	29/32	29/32	150
0100806042	140	45/40	29/32	150
0100806044	150	29/32	29/32	200
0100806046	165	45/40	29/32	200



General Information

Product name: FLAT FLANGE LIDS, 2 NECKS

Rev. Date: Jul 11, 2023

Details:

One central vertical neck, one side angled neck.

Picture



Technical Information

Ref.	Height (mm)	Center Neck	Side Neck	Laboratory Flange DN
0100806030	75	14/23	14/23	60
0100806032	85	29/32	14/23	60
0100806034	95	29/32	14/23	100

Product datasheet

Ref.	Height (mm)	Center Neck	Side Neck	Laboratory Flange DN
0100806036	95	29/32	29/32	100
0100806038	110	45/40	29/32	100
0100806040	125	29/32	29/32	150
0100806042	140	45/40	29/32	150
0100806044	150	29/32	29/32	200
0100806046	165	45/40	29/32	200

