

## Información general

**Nombre de producto:** TAPAS REACTORES 3 BOCAS

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

Una boca central vertical y 2 bocas laterales inclinadas.

## Fotografía



## Datos técnicos

Ref.	Altura (mm)	Boca Central	Bocas Laterales	DN Boca
0100806050	75	14/23	14/23	60
0100806052	85	29/32	14/23	60
0100806054	95	29/32	14/23	100

## Ficha de producto

Ref.	Altura (mm)	Boca Central	Bocas Laterales	DN Boca
0100806056	95	29/32	29/32	100
0100806058	110	45/40	29/32	100
0100806060	125	29/32	29/32	150
0100806062	140	45/40	29/32	150
0100806064	150	29/32	29/32	200
0100806066	165	45/40	29/32	200



## General Information

**Product name:** FLAT FLANGE LIDS, 3 NECKS

**Rev. Date:** Jul 11, 2023

**Details:**

One central vertical neck, two side angled necks.

## Picture



## Technical Information

Ref.	Height (mm)	Center Neck	Side Necks	Laboratory Flange DN
0100806050	75	14/23	14/23	60
0100806052	85	29/32	14/23	60
0100806054	95	29/32	14/23	100

# Product datasheet

Ref.	Height (mm)	Center Neck	Side Necks	Laboratory Flange DN
0100806056	95	29/32	29/32	100
0100806058	110	45/40	29/32	100
0100806060	125	29/32	29/32	150
0100806062	140	45/40	29/32	150
0100806064	150	29/32	29/32	200
0100806066	165	45/40	29/32	200

