

Información general

Nombre de producto: EJE AGITADOR DE ÁNCORA

Fecha Rev.: Jul 11, 2023

Detalles:

- Químicamente inerte.
- Fácil de limpiar.
- Evita la contaminación cruzada.
- Resistente a roturas y arañazos, y hasta 280 °C.
- Los modelos de eje expuesto (EX) tienen un núcleo de acero inoxidable, con un extremo expuesto de longitud de 50 mm para una fijación más rígida en el mandril.

Fotografía



Datos técnicos

Ref.	Ø Eje (mm)	Ø Rotor (mm)	Longitud (mm)
10AGAN0610	6	80	500
10AGAN0815	8	80	500
10AGAN1040	10	100	550
10AGAN1045	10	140	650EX
10AGAN1250	12	140	650EX
10AGAN1655	16	180	750EX
10AGAN1660	16	180	1000EX



General Information

Product name: SHAFT STIRRER, ANCHOR

Rev. Date: Jul 11, 2023

Details:

- Chemically inert.
- Easy to clean.
- No cross contamination.
- Resistant to breakage, scratches and heat up to 280 °C.
- Exposed stainless steel core models (EX) provide extra rigid clamping in the drive chuck, with an exposed length of 50 mm.

Picture



Technical Information

Ref.	Shaft Ø (mm)	Rotor Ø (mm)	Length (mm)
10AGAN0610	6	80	500
10AGAN0815	8	80	500
10AGAN1040	10	100	550
10AGAN1045	10	140	650EX
10AGAN1250	12	140	650EX
10AGAN1655	16	180	750EX
10AGAN1660	16	180	1000EX

