

Información general

Nombre de producto: GUÍAS PARA AGITADORES

Fecha Rev.: Jul 11, 2023

Detalles:

Compatible con todos los agitadores estándar.

Alta resistencia; permite trabajar con medio y alto vacío. Se adapta a bocas de vidrio normalizadas con depósito de refrigeración.

PTFE.

Fotografía



Datos técnicos

Ref.	Ø (MM)	Para esmerilado NS
10GUIA1926	6	19/26

Ficha de producto

Ref.	Ø (MM)	Para esmerilado NS
10GUIA1928	8	19/26
10GUIA2426	6	24/29
10GUIA2428	8	24/29
10GUIA2429	10	24/29
10GUIA2932	10	29/32
10GUIA2933	12	29/32
10GUIA3435	10	34/35
10GUIA3436	12	34/35
10GUIA3437	16	34/35
10GUIA4540	16	45/40



General Information

Product name: STIRRER GUIDES

Rev. Date: Jul 11, 2023

Details:

Compatible with all standard stirrers.

High resistance; suitable for medium and high vacuum. It can be used with normalized glass necks with cooling tank.

PTFE.

Picture



Technical Information

Ref.	Ø (MM)	Joint
10GUIA1926	6	19/26

Product datasheet

Ref.	∅ (MM)	Joint
10GUIA1928	8	19/26
10GUIA2426	6	24/29
10GUIA2428	8	24/29
10GUIA2429	10	24/29
10GUIA2932	10	29/32
10GUIA2933	12	29/32
10GUIA3435	10	34/35
10GUIA3436	12	34/35
10GUIA3437	16	34/35
10GUIA4540	16	45/40

