

Información general

Nombre de producto: EMBUDOS PESASUSTANCIAS

Fecha Rev.: Jul 11, 2023

Detalles:

Fabricados en vidrio borosilicato 3.3.

Ideal para el pesaje de muestras pequeñas.

Los reactivos en polvo pueden verterse o enjuagarse a través de la rama de los embudos de pesaje hacia el recipiente de recepción.

Fotografía



Datos técnicos

Ficha de producto

Ref.	Capacidad Aprox. (ML)	Para matraces aforados con tapón nº	Ø Exterior vástago (MM)
10EMPE0000	0.5	7/16	5
10EMPE0001	1	7/16	6
10EMPE0003	3	10/19	8
10EMPE0006	6	14/23	8
10EMPE0010	10	19/26	8
10EMPE0020	20	19/26	12



General Information

Product name: WEIGHING FUNNELS

Rev. Date: Jul 11, 2023

Details:

Made from borosilicate glass 3.3.

Ideal for weighing small samples.

Powder reagents can be poured or washed through the funnel stem towards a receiver container.

Picture



Technical Information

Ref.	Capacity Approx. (ML)	For volumetric flasks with stopper	OD Stem (MM)
10EMPE0000	0.5	7/16	5

Product datasheet

Ref.	Capacity Approx. (ML)	For volumetric flasks with stopper	OD Stem (MM)
10EMPE0001	1	7/16	6
10EMPE0003	3	10/19	8
10EMPE0006	6	14/23	8
10EMPE0010	10	19/26	8
10EMPE0020	20	19/26	12

