

Información general

Nombre de producto: CALCÍMETRO DE BERNARD

Fecha Rev.: Jul 11, 2023

Detalles:

Para la determinación porcentual del contenido de CaCO_3 (carbonato cálcico) según el método de Bernard.

La muestra se coloca en el matraz proporcionado, junto con el ácido clorhídrico vertido en un tubo de muestra. Tras reaccionar el ácido clorhídrico con la muestra, el dióxido de carbono formado por la reacción empuja el agua de la bureta graduada hacia la ampolla de nivel.

Fotografía



Datos técnicos

Ref.	Capacidad (ML)
10CABE6256	100

General Information

Product name: BERNARD CALCIMETER

Rev. Date: Jul 11, 2023

Details:

For the percentage determination of the CaCO_3 (calcium carbonate) content according to the De Bernard method.

The sample is placed in the flask provided, together with hydrochloric acid poured into a sample tube. After the hydrochloric acid has reacted with the sample, the carbon dioxide formed by the reaction pushes the water from the graduated burette into the glass level container.

Picture



Technical Information

Ref.	Capacity (ML)
10CABE6256	100 ML