

Información general

Nombre de producto: PAR GUANTES PROTECCIÓN

Fecha Rev.: Oct 17, 2024

Detalles:

Guantes de algodón suaves y flexibles. Excelentes para manipular objetos calientes recién salidos de la autoclave o del horno, o para otras protecciones contra el calor en el laboratorio, la escuela o la industria.

Marcado CE.

Resistentes al calor hasta 232°C (450°F); reemplazan los guantes peligrosos de asbesto.

Los guanteletes en dos longitudes brindan protección para la muñeca y el antebrazo.

El color naranja brillante ofrece alta visibilidad.

Nota: Los guantes no son impermeables y no deben usarse para manipular hielo seco.

Fotografía



Datos técnicos

| Ref. | Medidas (MM) |
|------------|-----------------------------------|
| 10GUAP0280 | Largo guante 180, largo total 470 |



General Information

Product name: PROTECTION GLOVES, PAIR

Rev. Date: Jul 11, 2023

Details:

Soft, pliant, all-cotton terry cloth gloves are superb for handling hot objects out of the autoclave or oven, or for other heat protection in the lab, school or industry.

CE marked.

Heat-resistant up to 232°C (450°F); replace hazardous asbestos gloves.

Gauntlets in two lengths provide wrist and forearm protection.

Bright orange color offers high visibility.

Note: Gloves are not waterproof and should not be used to handle dry ice.

Picture



LABOR

Technical Information

| Ref. | Dimensions (MM) |
|------------|--------------------------------------|
| 10GUAP0280 | Wrist length 280, mid-arm length 470 |

