

Información general

Nombre de producto: JUNTA SILICONA CON PTFE PARA TAPONES ROSCA GL

Fecha Rev.: Nov 05, 2024

Detalles:

Anillo de silicona con una arandela de PTFE. Una vez montado, el producto sólo queda expuesto al PTFE. Resistencia a la temperatura -60°C a 230°C.

Fotografía



Datos técnicos

Ref.	GL	Ø Orificio (mm)
0200162102	14	6
0200162108	18	10
0200162104	18	6
0200162106	18	8
0200162112	25	10
0200162114	25	12
0200162110	25	8
0200162116	32	10
0200162118	32	12
0200162120	32	14
0200162122	32	16
0200162124	32	18



General Information

Product name: SILICONE GASKET WITH PTFE FOR GL THREADED PLUGS

Rev. Date: Nov 05, 2024

Details:

Silicone ring with PTFE washer. After assembly, the product is only exposed to PTFE.
Temperature resistance from -60°C to 230°C.

Picture



Technical Information

Ref.	GL	Ø Hole (mm)
0200162102	14	6
0200162108	18	10
0200162104	18	6
0200162106	18	8
0200162112	25	10
0200162114	25	12
0200162110	25	8
0200162116	32	10
0200162118	32	12
0200162120	32	14
0200162122	32	16
0200162124	32	18

