

## Información general

**Nombre de producto:** GRADILLAS AGUJERO CUADRADO

**Fecha Rev.:** Nov 5, 2024

**Detalles:**

Apilable, ligera y ocupa poco espacio. Apropriada para su calentamiento en baño de agua, así como el almacenamiento de muestras en la nevera e incubación de estufas. Con posiciones marcadas con número y letra. Las gradillas se suministran desplegadas de forma plana y pueden armarse con pocas manipulaciones de forma firme e inseparable.

Temperatura de uso de -20 a 90 °C.

Dimensiones: 245 x 105 x 64 mm

## Fotografía



## Datos técnicos

Ref.	Tubo diámetro (MM)	Capacidad de tubos (Nº Tubos)
10GPCB2130	30	21
10GPCB2425	25	24
10GPCB4020	20	40
10GPCB6016	16	60
10GPCB9013	13	90

## General Information

**Product name:** TEST TUBE RACKS, SQUARE HOLES

**Rev. Date:** Nov 5, 2024


**Details:**

Stackable, lightweight and space-saving. Suitable for heating in a water bath, as well as storage of samples in the refrigerator and incubation ovens. With positions marked with number and letter. The racks are supplied flat unfolded and can be assembled with little handling firmly and inseparably.

Operating temperature -20 to 90 °C.

Dimensions: 245 x 105 x 64 mm

## Picture



## Technical Information

Ref.	Tube diameter (MM)	Tubes Capacity (N° Tubes)
10GPCB2130	30	21
10GPCB2425	25	24
10GPCB4020	20	40
10GPCB6016	16	60
10GPCB9013	13	90