

Información general

Nombre de producto: TAPÓN CON ORIFICIO A ROSCA RODAVISS SIN JUNTA

Fecha Rev.: Nov 19, 2024

Detalles:

Diseñado específicamente para trabajos que requieran condiciones de vacío y presión. Este tapón destaca por su capacidad de operar de manera segura a temperaturas de hasta 200 °C, lo que lo hace ideal para aplicaciones que implican altas temperaturas. Además, su sistema de unión de juntas esmeriladas y facilita la separación una vez finalizado el proceso químico minimiza el riesgo de rotura, aumentando la durabilidad y la seguridad durante su uso. Otra de sus ventajas es su capacidad de ser intercambiado con uniones esmeriladas conforme a la norma DIN 12 249.

Este producto debe ir acompañado de:

[Anillo de seguridad](#)

[Junta de nitrilo](#)

Fotografía



Datos técnicos

Ref.	Rosca
10TAOR0014	NS-14
10TAOR0019	NS-19
10TAOR0024	NS-24
10TAOR0029	NS-29
10TAOR0034	NS-34
10TAOR0045	NS-45



General Information

Product name: SCREW CAP WITH HOLE, RODAVISS, WITHOUT LINER

Rev. Date: Nov 15, 2024

Details:

Designed specifically for jobs requiring vacuum and pressure conditions. The plug is noted for its ability to operate safely at temperatures up to 200°C, making it ideal for applications involving high temperatures. In addition, its ground jointing system and easy separation after chemical processing minimises the risk of breakage, increasing durability and safety during use. Another advantage is its ability to be interchanged with ground joints in accordance with DIN 12 249.

This product must be accompanied by:

[Unlocking seal for rodaviss](#)

[Sealing nitrile ring](#)

Picture



Technical Information

Ref.	Thread
10TAOR0014	NS-14
10TAOR0019	NS-19
10TAOR0024	NS-24
10TAOR0029	NS-29
10TAOR0034	NS-34
10TAOR0045	NS-45

