


Información general

Nombre de producto: MATRACES FONDO REDONDO 3 BOCAS PARALELAS

Fecha Rev.: Dic 12, 2024

Detalles:

Fotografía



Datos técnicos

| Ref.       | Capacidad (ML) | Esmerilado P | Esmerilado N | Esmerilado O |
|------------|----------------|--------------|--------------|--------------|
| 0100109001 | 100            | 19/26        | 19/26        | 24/29        |
| 0100109002 | 100            | 14/23        | 14/23        | 29/32        |
| 0100109003 | 100            | 19/26        | 19/26        | 29/32        |
| 0100109004 | 100            | 29/32        | 29/32        | 29/32        |
| 0100109005 | 250            | 19/26        | 19/26        | 24/29        |
| 0100109006 | 250            | 14/23        | 14/23        | 29/32        |
| 0100109007 | 250            | 19/26        | 19/26        | 29/32        |
| 0100109008 | 250            | 29/32        | 29/32        | 29/32        |

Ficha de producto

| Ref.       | Capacidad (ML) | Esmerilado P | Esmerilado N | Esmerilado O |
|------------|----------------|--------------|--------------|--------------|
| 0100109009 | 500            | 19/26        | 19/26        | 24/29        |
| 0100109010 | 500            | 14/23        | 14/23        | 29/32        |
| 0100109011 | 500            | 19/26        | 19/26        | 29/32        |
| 0100109012 | 500            | 29/32        | 29/32        | 29/32        |
| 0100109013 | 1000           | 19/26        | 19/26        | 24/29        |
| 0100109014 | 1000           | 14/23        | 14/23        | 29/32        |
| 0100109015 | 1000           | 19/26        | 19/26        | 29/32        |
| 0100109016 | 1000           | 29/32        | 29/32        | 29/32        |
| 0100109017 | 1000           | 29/32        | 29/32        | 45/40        |
| 0100109018 | 2000           | 19/26        | 19/26        | 24/29        |
| 0100109019 | 2000           | 14/23        | 14/23        | 29/32        |
| 0100109020 | 2000           | 19/26        | 19/26        | 29/32        |
| 0100109021 | 2000           | 29/32        | 29/32        | 29/32        |
| 0100109022 | 2000           | 29/32        | 29/32        | 45/40        |
| 0100109023 | 3000           | 14/23        | 14/23        | 29/32        |
| 0100109024 | 3000           | 29/32        | 29/32        | 29/32        |
| 0100109025 | 3000           | 29/32        | 29/32        | 45/40        |
| 0100109026 | 4000           | 14/23        | 14/23        | 29/32        |
| 0100109027 | 4000           | 29/32        | 29/32        | 29/32        |
| 0100109028 | 4000           | 29/32        | 29/32        | 45/40        |
| 0100109029 | 5000           | 29/32        | 29/32        | 29/32        |
| 0100109030 | 5000           | 29/32        | 29/32        | 45/40        |
| 0100109031 | 6000           | 14/23        | 14/23        | 29/32        |
| 0100109032 | 6000           | 29/32        | 29/32        | 29/32        |
| 0100109033 | 6000           | 29/32        | 29/32        | 45/40        |



Ficha de producto

| Ref.       | Capacidad (ML) | Esmerilado P | Esmerilado N | Esmerilado O |
|------------|----------------|--------------|--------------|--------------|
| 0100109034 | 10000          | 29/32        | 29/32        | 29/32        |
| 0100109035 | 10000          | 29/32        | 29/32        | 45/40        |
| 0100109036 | 20000          | 29/32        | 29/32        | 45/40        |




General Information

Product name: MATRACES FONDO REDONDO 3 BOCAS PARALELAS

Rev. Date: Dec 12, 2024

Details:

Picture



Technical Information

| Ref.       | Capacity (ML) | Joint P | Joint N | Joint O |
|------------|---------------|---------|---------|---------|
| 0100109001 | 100           | 19/26   | 19/26   | 24/29   |
| 0100109002 | 100           | 14/23   | 14/23   | 29/32   |
| 0100109003 | 100           | 19/26   | 19/26   | 29/32   |
| 0100109004 | 100           | 29/32   | 29/32   | 29/32   |
| 0100109005 | 250           | 19/26   | 19/26   | 24/29   |
| 0100109006 | 250           | 14/23   | 14/23   | 29/32   |
| 0100109007 | 250           | 19/26   | 19/26   | 29/32   |
| 0100109008 | 250           | 29/32   | 29/32   | 29/32   |

## Product datasheet

| Ref.       | Capacity (ML) | Joint P | Joint N | Joint O |
|------------|---------------|---------|---------|---------|
| 0100109009 | 500           | 19/26   | 19/26   | 24/29   |
| 0100109010 | 500           | 14/23   | 14/23   | 29/32   |
| 0100109011 | 500           | 19/26   | 19/26   | 29/32   |
| 0100109012 | 500           | 29/32   | 29/32   | 29/32   |
| 0100109013 | 1000          | 19/26   | 19/26   | 24/29   |
| 0100109014 | 1000          | 14/23   | 14/23   | 29/32   |
| 0100109015 | 1000          | 19/26   | 19/26   | 29/32   |
| 0100109016 | 1000          | 29/32   | 29/32   | 29/32   |
| 0100109017 | 1000          | 29/32   | 29/32   | 45/40   |
| 0100109018 | 2000          | 19/26   | 19/26   | 24/29   |
| 0100109019 | 2000          | 14/23   | 14/23   | 29/32   |
| 0100109020 | 2000          | 19/26   | 19/26   | 29/32   |
| 0100109021 | 2000          | 29/32   | 29/32   | 29/32   |
| 0100109022 | 2000          | 29/32   | 29/32   | 45/40   |
| 0100109023 | 3000          | 14/23   | 14/23   | 29/32   |
| 0100109024 | 3000          | 29/32   | 29/32   | 29/32   |
| 0100109025 | 3000          | 29/32   | 29/32   | 45/40   |
| 0100109026 | 4000          | 14/23   | 14/23   | 29/32   |
| 0100109027 | 4000          | 29/32   | 29/32   | 29/32   |
| 0100109028 | 4000          | 29/32   | 29/32   | 45/40   |
| 0100109029 | 5000          | 29/32   | 29/32   | 29/32   |
| 0100109030 | 5000          | 29/32   | 29/32   | 45/40   |
| 0100109031 | 6000          | 14/23   | 14/23   | 29/32   |
| 0100109032 | 6000          | 29/32   | 29/32   | 29/32   |
| 0100109033 | 6000          | 29/32   | 29/32   | 45/40   |



Product datasheet

| Ref.       | Capacity (ML) | Joint P | Joint N | Joint O |
|------------|---------------|---------|---------|---------|
| 0100109034 | 10000         | 29/32   | 29/32   | 29/32   |
| 0100109035 | 10000         | 29/32   | 29/32   | 45/40   |
| 0100109036 | 20000         | 29/32   | 29/32   | 45/40   |

