

Información general

Nombre de producto: BURETAS CON DEPÓSITO CLASE A

Fecha Rev.: Jul 11, 2023

Detalles:

Con franja azul sobre fondo blanco.

Fabricado en vidrio borosilicato 3.3.

Fotografía



Datos técnicos

Ref.	Llave	Capacidad depósito (ML)	Capacidad (ML)	Subdivisiones (ML)
0100490602	Vidrio	250	2	0.01
0100490603	Vidrio	250	5	0.02

Ficha de producto

Ref.	Llave	Capacidad depósito (ML)	Capacidad (ML)	Subdivisiones (ML)
0100490604	Vidrio	250	10	0.05
0100490607	Vidrio	250	25	0.1
0100490608	Vidrio	500	50	0.1
0100490612	PTFE Cónica	250	2	0.01
0100490613	PTFE Cónica	250	5	0.02
0100490614	PTFE Cónica	250	10	0.05
0100490617	PTFE Cónica	250	25	0.1
0100490618	PTFE Cónica	500	50	0.1
0100490622	PTFE Aguja	250	2	0.01
0100490623	PTFE Aguja	250	5	0.02
0100490624	PTFE Aguja	250	10	0.05
0100490627	PTFE Aguja	250	25	0.1
0100490628	PTFE Aguja	500	50	0.1

General Information

Product name: BURETTES WITH RESERVOIR, CLASS A

Rev. Date: Jul 11, 2023

Details:

With blue strip on white background.

Made from borosilicate glass 3.3.

Picture



Technical Information

Ref.	Stopcock	Reservoir Capacity (ML)	Capacity (ML)	Subdivisiones (ML)
0100490602	GLASS	250	2	0.01
0100490603	GLASS	250	5	0.02

Product datasheet

Ref.	Stopcock	Reservoir Capacity (ML)	Capacity (ML)	Subdivisiones (ML)
0100490604	GLASS	250	10 ML	0.05
0100490607	GLASS	250	25 ML	0.1
0100490608	GLASS	500	50 ML	0.1
0100490612	CONICAL PTFE	250	2	0.01
0100490613	CONICAL PTFE	250	5 ML	0.02
0100490614	CONICAL PTFE	250	10 ML	0.05
0100490617	CONICAL PTFE	250	25 ML	0.1
0100490618	CONICAL PTFE	500	50 ML	0.1
0100490622	PTFE NEEDLE	250	2	0.01
0100490623	PTFE NEEDLE	250	5 ML	0.02
0100490624	PTFE NEEDLE	250	10 ML	0.05
0100490627	PTFE NEEDLE	250	25 ML	0.1
0100490628	PTFE NEEDLE	500	50 ML	0.1