

## Información general

**Nombre de producto:** DISPENSADOR DE LÍQUIDOS

**Fecha Rev.:** Jul 11, 2023

**Detalles:**

Para dispensar líquidos con precisión y repetitividad.

## Fotografía



## Datos técnicos

Ref.	Rango Volumen (mL)	Tolerancia ( $\pm$ ML)	Subdivisiones (ML)
10DOFS0002	0.25-2.5	$\pm 0.015$	0.05
10DOFS0005	0.5-5	$\pm 0.03$	0.1
10DOFS0010	1-10	$\pm 0.060$	0.2

## Ficha de producto

Ref.	Rango Volumen (mL)	Tolerancia ( $\pm$ ML)	Subdivisiones (ML)
10DOFS0030	2.5-30	$\pm 0.18$	0.5
10DOFS0060	5-60	$\pm 0.36$	1
10DOFS0100	10-100	$\pm 0.6$	2



## General Information

**Product name:** BOTTLE TOP DISPENSER

**Rev. Date:** Jul 11, 2023

**Details:**

Precise and repeatable liquid dispensing.

## Picture



## Technical Information

Ref.	Volume Range (mL)	Tolerance ( $\pm$ ML)	Subdivisiones (ML)
10DOFS0002	0.25-2.5	$\pm$ 0.015	0.05
10DOFS0005	0.5-5	$\pm$ 0.03	0.1
10DOFS0010	1-10	0.06 ML	0.2

## Product datasheet

Ref.	Volume Range (mL)	Tolerance ( $\pm$ ML)	Subdivisiones (ML)
10DOFS0030	2.5-30	$\pm 0.18$	0.5
10DOFS0060	5-60	$\pm 0.36$	1
10DOFS0100	10-100	$\pm 0.6$	2

