

Información general

Nombre de producto: MICROPIPETAS MONOCANAL DE VOLUMEN VARIABLE - HTL MODELO DISCOVERY COMFORT

Fecha Rev.: Jul 11, 2023

Detalles:

Combinan el mayor nivel de precisión con un diseño tradicional para que el pipeteado sea mas cómodo y seguro.

Correcto agarre reduciendo el estrés en la muñeca y la fatiga, causantes de lesión por el esfuerzo repetitivo.

Diseño robusto.

Códigos de color

Resistente a luz UV y totalmente autoclavables.

Expulsor de punta secuencial y con mínimo esfuerzo.

Fotografía



Datos técnicos

Ref.	Rango Volumen (μL)	Color
154050D250	50-250	Verde
154051D2	0.2-2.0	Rojo
154052D10	0,5-10	Rojo
154053D20	2-20	Amarillo
154054D100	10-100	Amarillo
154055D200	20-200	Amarillo
154056D1000	100-1000	Azul
154057D50	5-50	Amarillo
154058D5000	500-5000	Blanco
154059D10000	1000-10000	Blanco



General Information

Product name: HTL SINGLE CHANNEL ADJUSTABLE VOLUME MICROPIPETTE, MODEL DISCOVERY COMFORT

Rev. Date: Jul 11, 2023

Details:

Combination of the highest level of accuracy and precision with traditional design for a more comfortable and safe pipetting.

Its ergonomic grip reduces stress on the wrist and fatigue, thus preventing repetitive strain injuries.

Robust design.

Colour coding.

UV resistant and fully autoclavable.

Sequential tip ejection with minimal effort.

Picture



Technical Information

Ref.	Range Volume (μ L)	Colour
154050D250	50-250	Verde
154051D2	0.2-2.0	Red
154052D10	0.5-10	Red
154053D20	2-20	Amarillo
154054D100	10-100	Amarillo
154055D200	20-200	Amarillo
154056D1000	100-1000	Blue
154057D50	5-50	Amarillo
154058D5000	500-5000	White
154059D10000	1000-10000	White

