

Información general

Nombre de producto: MICROPIPETAS MULTICANAL DE VOLUMEN VARIABLE - HTL MODELO DISCOVERY COMFORT

Fecha Rev.: Jul 11, 2023

Detalles:

Disponible en 8 y 12 canales.

Excelente precisión en series de pipeteado de larga duración.

Diseño robusto.

Sellado de la punta seguro y sencillo.

Código de color.

Resistente a luz UV y totalmente autoclavable.

Expulsor de punta secuencial y con mínimo esfuerzo.

Fotografía



BR
3
0
1

Datos técnicos

Ref.	Rango Volumen (μL)	Color	Canales
155121DV8-10	1-10	Rojo	8
155122DV8-50	5-50	Amarillo	8
155123DV8-200	20-200	Amarillo	8
155124DV8-300	50-300	Verde	8
155125DV12-10	1-10	Rojo	12
155126DV12-50	5-50	Amarillo	12
155127DV12200	20-200	Amarillo	12
155128DV12300	50-300	Verde	12



General Information

Product name: HTL MULTICHANNEL ADJUSTABLE VOLUME MICROPIPETTE, MODEL DISCOVERY COMFORT

Rev. Date: Jul 11, 2023

Details:

Available in 8 and 12 channel version.

Excellent precision in long pipetting series.

Robust design.

Secure and easy tip sealing.

Colour coding.

UV resistant and fully autoclavable.

Sequential tip ejection with minimal effort.

Picture



Technical Information

Ref.	Range Volume (μ L)	Colour	Channels
155121DV8-10	1-10	Red	8
155122DV8-50	5-50	Amarillo	8
155123DV8-200	20-200	Amarillo	8
155124DV8-300	50-300	Verde	8
155125DV12-10	1-10	Red	12
155126DV12-50	5-50	Amarillo	12
155127DV12200	20-200	Amarillo	12
155128DV12300	50-300	Verde	12

