

Información general

Nombre de producto: MICROPIPETAS DE VOLUMEN VARIABLE AUTOCLAVABLE LABGILLS

Fecha Rev.: Jul 11, 2023

Detalles:

Totalmente autoclavables.

Diseño ergonómico y peso ligero.

Selección de volumen mediante contador con parada con clic.

Calibración de acuerdo a las normas DIN 12650 e ISO 8655, asegurando una elevada exactitud y precisión.

Fotografía



Datos técnicos

Ficha de producto

Ref.	Rango Volumen (μ L)	Incremento (μ l)
18MVLG00002	0,1-2,5	0,01
18MVLG00010	0,5-10	0,1
18MVLG00020	2-20	0,1
18MVLG00050	5-50	0,5
18MVLG00100	10-100	0,5
18MVLG00200	20-200	1,0
18MVLG01000	100-1000	5,0
18MVLG05000	500-5000	50,0
18MVLG10000	500-10000	100,0



General Information

Product name: LABGILLS ADJUSTABLE VOLUME MICROPIPETTES, AUTOCLAVABLE

Rev. Date: Jul 11, 2023

Details:

Totally autoclavable.

Ergonomic design and lightweight.

Volume selection through click counter.

Calibration according to DIN 12650 and ISO 8655 standards, ensuring high precision and accuracy.

Picture



Technical Information



Product datasheet

Ref.	Range Volume (μL)	Increment (μL)
18MVLG00002	0.1-2.5	0.01
18MVLG00010	0.5-10	0.1
18MVLG00020	2-20	0.1
18MVLG00050	5-50	0.5
18MVLG00100	10-100	0.5
18MVLG00200	20-200	1.0
18MVLG01000	100-1000	5.0
18MVLG05000	500-5000	50.0
18MVLG10000	500-10000	100.0

