

Información general

Nombre de producto: ADAPTADOR M/M ACODADO A 75°

Fecha Rev.: Jul 11, 2023

Detalles:

Fotografía



Datos técnicos

Ref.	Esmerilado Macho
0100220111	14/23
0100220122	19/26
0100220133	24/29
0100220144	29/32

General Information

Product name: ADAPTER 2 CONES, BENT 75°

Rev. Date: Jul 11, 2023

Details:

Picture



Technical Information

Ref.	Cone Joint
0100220111	14/23
0100220122	19/26
0100220133	24/29
0100220144	29/32